


BVI

A Bright Vision on Illumination

Light Fibers




Standard
Light Fibers



Wide View
Light Fibers



45 Degree Shielded
Light Fibers



Vitreco

Vitreoretinal Surgical Products

bvimedical.com

Light Fibers

Light Fibers featuring high-end specifications with consistent quality performance for ophthalmic endo-illumination procedures.

Rigid shaft

- New hardened metal shaft reduces flexibility by 25%
- Reduced bend ability improves surgical control

Ergonomically designed

- Ergonomic handle design for improved intraocular surgical control and precise surgical manoeuvres
- Compatible with most endo-illumination systems due to light source adaptor




Light fiber tip design

- In Standard, Wide view and 45 Degree Shielded tip designs


Wide range of sizes

- Available in 23G, 25G and 27G with colour coding for identifying gauge size instrument




Standard Light Fibers

Straight tip for focal endo-illumination.



Wide View Light Fibers

Wide view tip for global endo-illumination.



45 Degree Shielded Light Fibers

45 Degree shielded wide view tip for shielded endo-illumination.

Light Fibers Options



Standard Light Fibers

Description	SKU #
● 23G Standard light fibers	LF23.D01
● 25G Standard light fibers	LF25.D01
● 27G Standard light fibers	LF27.D01



45 Degree Shielded Light Fibers

Description	SKU #
● 23G 45 Degree shielded light fibers	LF23.D03
● 25G 45 Degree shielded light fibers	LF25.D03
● 27G 45 Degree shielded light fibers	LF27.D03



Wide View Light Fibers

Description	SKU #
● 23G Wide view light fibers	LF23.D02
● 25G Wide view light fibers	LF25.D02
● 27G Wide view light fibers	LF27.D02

Gauge Size with Colour Coding



Packaged box/5, Sterile



Light Source Adaptors

The Light Fibers are compatible with most endo-illumination systems in combination with Light Source Adaptors*.

Packaged box/1, Sterile

Vitreq® Light Source Adaptor	Description	SKU #
	Constellation® Light Source Adaptor	LF00.R01
	Stellaris® Light Source Adaptor	LF00.R02
	Geuder® Light Source Adaptor	LF00.R03
	Ruck® Light Source Adaptor	LF00.R04
	DORC® Light Source Adaptor	LF00.R05
	Vitreq® Light Fiber connector	LF23.D01

* All product names, trademarks and registered trademarks are property of their respective owners.

FOR COMPLETE PRODUCT INSTRUCTIONS, PLEASE REFER TO THE BVI IFU OR PRODUCT INSERT

Questions: bvimedical.com/customer-support

