

BVI

Fluid Control



PVA Eye Spears

Spear overall handle length is 7.16 cm or 2.82 in.

UltraClean™ LASIK PVA Eye Spears

Compressed shape, designed for fast expanding. Made from highly absorbent, virtually lint free PVA material. Designed specifically for rapid recovery of fluids during LASIK and other corneal procedures.



Qty: 100/Bx
Brand: UltraClean™

Description	SKU #
UltraClean™ Eye Spears	40810
UltraClean™ Eye Spears	40815

UltraClean™ LASIK PVA Expanded Eye Spears

Designed for instantaneous wicking, pre-expanded tip. Made from highly absorbent, virtually lint free PVA material. No need to pre-soak and fan.



Qty: Various/Bx
Brand: UltraClean™

Description	SKU #
UltraClean™ PVA Expanded LASIK Eye Spear, 40/Bx	40800
UltraClean™ LASIK Expanded Eye Spear, 100/bx	40805

Visitec® LASIK PVA Eye Spears

Made from highly absorbent PVA material. Designed for optimal tissue manipulation and controlled fluid absorption to facilitate flap alignment and adhesion.



Qty: 100/Bx
Brand: Visitec®

Description	SKU #
LASIK Eye Spear - Orange	581761

Visitec® LASIK PVA Expanded Eye Spears

Made from highly absorbent PVA material. Pre-expanded for convenience. Designed to provide instantaneous wicking. No need to pre-soak and fan.



Qty: 40/Bx
Brand: Visitec®

Description	SKU #
LASIK Expanded Spear	581709

PVA Eye Spears

Made from exclusive polyvinyl acetal (PVA) material. The material is designed to be highly absorbent, fast wicking and virtually lint free.



Qty: Various/Bx
Brand: BVI®

Description	SKU #
PVA Eye Spear, 10 Spears/Pouch, 180/Bx	400101
PVA Eye Spear, 5 Spears/Pouch, 50/Bx	405101

Ultracell® PVA Eye Spears

Made from highly absorbent PVA material. Designed for rapid recovery of fluids during cataract and other general ophthalmic procedures.



Qty: Various/Bx
Brand: Ultracell®

Description	SKU #
PVA Eye Spear with dark blue handle, 100/Bx	40400-8
PVA Eye Spear with dark blue handle, 200/Bx	40405-8

Cellulose Eye Spears

Ultracell® Cellulose Eye Spears

Made from highly absorbent, natural, cellulose material. Fast wicking action delivers quick absorption of blood and fluid from the surgical area. The flat tip provides precise application at the surgical incision.



Qty: Various/Bx
Brand: Ultracell®

Description	SKU #
Ultracell® Cellulose Eye Spear with light blue handle, 144/Bx	40413
Ultracell® Cellulose Eye Spear with light blue handle, 200/Bx	40415

Visitec® Cellulose Eye Spears

Made from highly absorbent, natural, cellulose material. Fast wicking action delivers quick absorption of fluid from the surgical area. The flat tip provides precise application at the surgical incision. Magenta handle.



Qty: 100/Bx
Brand: Visitec®

Description	SKU #
Sponge Eye, 7 cm (2.75 in) overall length.	581089

Weck-Cel® Cellulose Eye Spears

Made from highly absorbent, natural, cellulose material. Designed for use in delicate, surgical areas and to maintain rigidity during wicking process. Light blue handle.



Qty: Various/Bx
Brand: Weck-Cel®

Description	SKU #
Weck-Cel® Spears, 6 Pack, 144/Bx	0008680
Weck-Cel® Spears, 10 Pack, 180/Bx	0008685

Sponge Points, Strips and Wicks

Weck-Cel® Sponge Points and Strips

Made from highly absorbent, natural, cellulose material.



Qty: Various/Bx
Brand: Weck-Cel®

Description	SKU #
Weck-Cel® Points, 240/Bx	0008661
Weck-Cel® Points and Strips, 24/Bx	0008660

Sponge Points and Eye Wicks

Made from exclusive polyvinyl acetal (PVA) material. The material is highly absorbent, fast wicking and virtually lint free.



Qty: Various/Bx
Brand: BVI®

Description	SKU #
Eye Wicks, 8 in, 2/Pouch, 20/Bx	400111
Sponge Points, 10 Points/Pouch, 200/Bx	410115
Sponge Points, 20/Pouch, 400/Bx	400115

Visitec® Eye Fluid Wick

Wicks absorb and direct fluid from surgical area into fluid collection bag.



Qty: 20/Bx
Brand: Visitec®

Description	SKU #
Eye fluid wick, 15 cm Long x 0.4 cm Wide	581048

Corneal Light Shields

Ultracell® Corneal Light Shield

Used to protect hinge and posterior side of corneal flap and protect the retina from the operating microscope light. Made from PVA material.




Qty: 20/Bx
Brand: Ultracell®

Description	SKU #
Ultracell® Corneal Light Shield, 7 mm (0.28 in)	40420

Keracel® Corneal Light Shield

Round shield cut flat at one end to fit against the hinge of a LASIK corneal flap. Made from a patented high-density formulation of virtually lint free Polyvinyl Acetal (PVA) sponge. 10 mm diameter with straight edge.



Qty: 20/Bx
Brand: Keracel®

Description	SKU #
Keracel® Corneal Light Shield with straight edge	400121

Corneal Light Shield

Used to protect hinge and posterior side of flap. Used to protect the retina from operating microscope light. The exclusive material is highly absorbent, fast wicking, virtually lint free, polyvinyl acetal (PVA).



Qty: 20/Bx
Brand: BVI®

Description	SKU #
Corneal Light Shield, 8 mm (0.3 in) diameter	400106

Visitec® Corneal Light Shield

PVA absorbent pad. Soft, porous structure designed to protect the retina from operating microscope light and contains no cellulose. Upon moistening, shield becomes soft and conforms to shape of cornea. Designed to keep cornea moist.



Qty: 20/Bx
Brand: Visitec®

Description	SKU #
Visitec® Corneal Light Shield, 8 mm (0.3 in) diameter	581062

Full Body Drapes


Visidrape® full body drapes feature wider widths and a three-sided channel bag for maximum fluid capacity. Full body drapes can be utilized in all ophthalmic procedures.

Full Body Incise Drape

Qty: 10/Bx

Brand: Visidrape®

Description	SKU #
Visidrape® Non-Woven Incise Drape with Fluid Collection Bag	581087
Overall size:	140 cm x 244 cm (55 in x 96 in)
Incise area:	9 cm x 11.5 cm (3 1/2 in x 4 1/2 in)
Fluid collection bag:	25.4 cm x 15.2 cm (10 in x 5.98 in)
Capacity:	500 ml
Material Type:	Non-Woven Wetlaid




Full Body Aperture Drape

Qty: 10/Bx

Brand: Visidrape®

Description	SKU #
Visidrape® Full Body Aperture Drape	581088
Overall size:	140 cm x 244 cm (55 in x 96 in)
Aperture area:	5 cm x 7.5 cm (2 in x 3 in)
Capacity:	500 ml
Material Type:	Non-Woven Wetlaid




Full Body Incise Drape with Plastic Head and Channel Bag

Qty: 10/Bx

Brand: Visidrape®

Description	SKU #
Visidrape® Full Body Incise Drape with Plastic Head and Channel Bag	581715
Overall size:	165 cm x 254 cm (65 in x 100 in)
Incise Area:	9 cm x 11 cm (3.5 in x 4 in)
Channel Bag:	62 cm x 40 cm (24 in x 15.75 in)
Capacity:	1000 ml
Material Type:	SMS (63 in) and Opaque (37 in)




Full Body Drapes

Non-Woven Full Body Aperture Drape

Qty: 10/Bx
Brand: Visidrape®

Description		SKU #
Visidrape® Non-Woven Full Body Aperture Drape		581722
Overall size:	165 cm x 254 cm (65 in x 100 in)	
Aperture area:	6 cm x 4 cm (2.25 in x 1.5 in)	
Material Type:	Non-Woven Wetlaid	



Half Body Drapes

Embossed Plastic Incise Drape with Fluid Collection Bag

Qty: 10/Bx
Brand: Visidrape®


Description		SKU #
Visidrape® Embossed Plastic Incise Drape with Fluid Collection Bag		581071
Overall size:	101.5 cm x 127 cm (40 in x 50 in)	
Incise area:	9 cm x 11.5 cm (3 1/2 in x 4 1/2 in)	
Fluid collection bag:	15.2 cm x 16.5 cm (6 in x 6 1/2 in)	
Capacity:	250 ml	
Material Type:	Polyethylene Opaque Film	



Embossed Plastic Aperture Drape with Fluid Collection Bag

Qty: 10/Bx
Brand: Visidrape®

Description		SKU #
Visidrape® Embossed Plastic Aperture Drape with Fluid Collection Bag		581078
Overall size:	101.5 cm x 127 cm (40 in x 50 in)	
Aperture area:	5 cm x 7.5 cm (2 in x 3 in)	
Fluid collection bag:	15.2 cm x 16.5 cm (6 in x 6 1/2 in)	
Capacity:	250 ml	
Material Type:	Polyethylene Opaque Film	




Half Body Drapes

Non-Woven Incise Drape with Fluid Collection Bag

Qty: 10/Bx
Brand: Visidrape®


Description		SKU #
Visidrape® Non-Woven Incise Drape with Fluid Collection Bag		581072
Overall size:	101.5 cm x 127 cm (40 in x 50 in)	
Incise area:	9 cm x 11.5 cm (3 1/2 in x 4 1/2 in)	
Fluid collection bag:	15.2 cm x 16.5 cm (6 in x 6 1/2 in)	
Capacity:	250 ml	
Material Type:	Non-Woven Wetlaid	



Non-Woven Aperture Drape with Fluid Collection Bag

Qty: 10/Bx
Brand: Visidrape®


Description		SKU #
Visidrape® Non-Woven Aperture Drape with Fluid Collection Bag		581079
Overall size:	101.5 cm x 127 cm (40 in x 50 in)	
Aperture area:	5 cm x 7.5 cm (2 in x 3 in)	
Fluid collection bag:	15.2 cm x 16.5 cm (6 in x 6 1/2 in)	
Capacity:	250 ml	
Material Type:	Non-Woven Wetlaid	



Non-Woven Incise Drape with Channel Bag

Qty: 10/Bx
Brand: Visidrape®

Description		SKU #
Visidrape® Non-Woven Incise Drape with Channel Bag		581084
Overall size:	140 cm x 127 cm (55 in x 50 in)	
Incise area:	9 cm x 11.5 cm (3 1/2 in x 4 1/2 in)	
Fluid collection bag with channels:	14.6 cm x 25.4 cm (5 3/4 in x 10 in)	
Capacity:	500 ml	
Material Type:	Non-Woven Wetlaid	




Half Body Drapes

Non-Woven Aperture Drape with Channel Bag

Qty: 10/Bx
Brand: Visidrape®


Description	SKU #
Visidrape® Non-Woven Aperture Drape with Channel Bag	581086
Overall size:	140 cm x 127 cm (55 in x 50 in)
Aperture area:	5 cm x 7.5 cm (2 in x 3 in)
Fluid collection bag with channels:	14.6 cm x 25.4 cm (5 3/4 in x 10 in)
Capacity:	500 ml
Material Type:	Non-Woven Wetlaid



Non-Woven Half Body Incise Drape* with Channel Bag

Qty: 10/Bx
Brand: Visidrape®

Description	SKU #
Visidrape® Non-Woven Half Body Incise Drape* with Channel Bag	581723
Overall size:	165 cm x 183 cm (65 in x 72 in)
Incise area:	9 cm x 11 cm (3.5 in x 4 in)
Channel bag:	62 cm x 40 cm (24 in x 15.75 in)
Capacity:	1000 ml
Material Type:	SMS



*Without wick

Mini Drapes

Ideal for refractive or office procedures.

Mini Plastic Incise Drape

Adhesive backing.

Qty: 10/Bx
Brand: Visidrape®

Description	SKU #
Visidrape® Mini Plastic Incise Drape	581095
Overall size:	41 cm x 41 cm (16 in x 16 in)
Incise area:	5 cm x 7.5 cm (2 in x 3 in)
Material Type:	Polyethylene Opaque Film



Mini Drapes

Mini Plastic Aperture Drape

Adhesive backing.

Qty: 10/Bx

Brand: Visidrape®

Description	SKU #
Visidrape® Mini Plastic Aperture Drape	581096
Overall size:	41 cm x 41 cm (16 in x 16 in)
Aperture area:	5 cm x 7.5 cm (2 in x 3 in)
Material Type:	Polyethylene Opaque Film




Mini Plastic Aperture Drape

Small incise provides complete coverage of lashes, lids and brows. Converts aperture to incise drape.

Qty: 10/Bx

Brand: Visidrape®

Description	SKU #
Visidrape® Incise Drape	581700
Overall size:	10.2 x 20.4 cm (4 x 8 in)
Adhesive area:	11 cm x 16 cm (4 in x 6 1/4 in)
Material Type:	Polyethylene Opaque Film



Dual Refractive Drape

Qty: 10/Bx

Brand: Visidrape®

Description	SKU #
Visidrape® Dual Refractive Drape	581710
Overall size:	40 cm x 26 cm (16 in x 10 1/4 in)
Aperture area:	2.5 cm x 3.75 cm (1 in x 1 1/2 in)
Material Type:	Polyethylene Opaque Film



Eye Drains

Keracel® Chayet Drainage Ring

Oval-shaped sponge to wick fluid and serve as a soft, virtually lint-free bed to place the corneal flap.



Qty: 20/Bx
Brand: Keracel®

Description	SKU #
Keracel® Chayet Drainage Ring, 12 mm	400122

Visitec® LASIK Eye Drain (Chayet)

Designed to reduce stromal bed hydration for a consistent ablation and for nasally hinged flaps.



Qty: 20/Bx
Brand: Visitec®

Description	SKU #
LASIK Eye Drain, 10 mm	581098

Fluid Collection

Eye Drain

Provides continuous control of fluid throughout the procedure. The wicks are made from highly absorbent, fast wicking, virtually lint free PVA material.



400108

Qty: Various/Bx
Brand: BVI®

Description	SKU #
Eye Drain with 80 cc Bag Bag: 21.6 cm x 5.1 cm (8.5 in x 2 in) Wick: Includes one, 15.2 cm (6 in)	400108
Eye Drain with 400 cc Bag Bag: 21.9 cm x 10.2 cm (8.6 in x 4 in) Wick: Includes two, 20.3 cm (8 in)	400109

Visitec® Fluid Collection Bags

For procedures where collection and control of fluid is designed to keep surgeon, patient and floor dry. Reduces operating room cleanup and turnaround time. Bags have an adhesive backing to adhere to surgical drape.



581005



581070

Qty: 20/Bx
Brand: Visitec®

Description	SKU #
Fluid Collection Bag: 17 cm x 16 cm (6 3/4 x 6 1/4 in), 250 ml capacity	581005
Fluid Collection Bag with Wick Bag: 17 cm x 16 cm (6 3/4 x 6 1/4 in), 250 ml capacity Wick: Includes one, 0.4 cm x 15 cm (5/32 x 5 7/8 in)	581006
Mini Fluid Collection with Wick Bag: 9 x 16 cm (3 1/2 x 6 1/4 in), 100 ml capacity Wick: Includes one, 0.4 x 15 cm (5/32 x 5 7/8 in)	581070

Instrument Wipes

Instrument Wipes

Wipes are designed to provide the surgeon a wipe for safe cleaning of fine microsurgical instruments. The exclusive material is highly absorbent, fast wicking, virtually lint free, polyvinyl acetal (PVA).



Qty: Various/Bx
Brand: BVI®

Description	SKU #
Instrument Wipe 92.1 x 76.2 mm (3.625 x 3 in), 20/Bx	400200
Wipe with Wick 76.2 x 76.2 mm (3 x 3 in), 20/Bx	223621
Wipe with Adhesive Back 50.8 x 50.8 mm (2 x 2 in), 20/Bx	400204

Hardten Collector

Adhesive-backed 50.8 mm x 50.8 mm (2 in x 2 in) sponge. May be attached to skin surface temporarily for fluid absorption. Sponge will hold up to 20 cc. Single-use, packaged sterile.



Qty: 20/Bx
Brand: BVI®

Description	SKU #
Hardten Collector	400126

Diamond Cleaning Block

4.5 cm block. Made of highly absorbent, virtually lint-free PVA material. Single-use, sterile.



Qty: 1/Bx
Brand: BVI®

Description	SKU #
Diamond Cleaning Block	400119

Visitec® Instrument Wipes



Qty: 20/Bx
Brand: Visitec®

Description	SKU #
Instrument Wipe 76.2 x 76.2 mm (3 x 3 in), 2 mm thick	90900

Visitec® Visiwipe Instrument Wipes



Qty: 20/Bx
Brand: Visitec®

Description	SKU #
Instrument Wipe 7.3 x 7.3 cm (2 7/8 x 2 7/8 in)	581047
Adhesive Wipe 7.6 x 4.5 cm (3 x 1 3/4 in)	581076

Other Fluid Control

Absorbent Stick

6.6 cm overall length x 5 mm diameter (2.59 in x 3/16 in).



Qty: 250/Bx
Brand: Visitec®

Description	SKU #
Absorbent Stick 6.6 cm, 10 per pouch, 25 pouches per box.	581060

Visiswab Cotton Swabs

7.6 cm (3 in) overall length. Rounded end is made for scleral depression and pointed end may be used as a second instrument to move sutures, manipulate tissue or absorb blood and fluid. Tightly wound ends create less linting than other cotton swabs. 10 per pouch, 25 pouches per box.



Qty: 250/Bx
Brand: Visitec®

Description	SKU #
Visiswab Cotton 7.6 cm (3 in) length	581069

Haddad Spatula

Spatula-shaped sponge that is malleable to any desired configuration. May be used to dry corneal bed, manipulate corneal tissue, or dry corneal surface trough. Single use, packaged sterile.




Qty: 20/Bx
Brand: BVI®

Description	SKU #
Haddad Spatula	400124

Visitec® Fluid Administration Set

Roller valve allows control of fluid flow. Sets with self-priming filter protect against particulates and microbiological contamination of 0.22 micron and larger.



Qty: 10/Bx
Brand: Visitec®

Description	SKU #
Set Single Fluid Administration Tubing length: 203 cm (80 in)	581014

Fluid Control Sets

Cataract Set



Qty: 10/Bx
Brand: BVI®

Description	SKU #
Cataract Kit	400113
Contents	SKU #
Eye spears, Qty/5	400101
Instrument wipe, Qty/1	400200
Corneal light shield, Qty/1	400106
80 cc Collection bag with wick, Qty/1	400108

Ophthalmic Surgery Set



Qty: 20/Bx
Brand: Visitec®

Description	SKU #
Ophthalmic Surgery Kit	581066
Contents	SKU #
Instrument wipe, Qty/1	581047
Eye fluid wick, Qty/1	581048
Corneal light shield, Qty/1	581062

Instrument Wipe Set



Qty: 20/Bx
Brand: BVI®

Description	SKU #
Instrument Wipe Set	223622
Contents	SKU #
Instrument wipe with wick, Qty/1	223621
Corneal light shield, Qty/1	400106

Non-Ophthalmic

Nasal Packs

Ultracell® brand PVA sponges for postoperative packing following nasal surgery.

- Used following septoplasty, rhinoplasty, and turbinectomy procedures
- May aid in holding septum firmly in place
- Available in a range of sizes: 3 cm pediatric to 10 cm
- Select sizes available in an anatomical shape
- Popular sizes also available in our patented layered design
- Excellent alternative to cumbersome gauze packing
- Latex-free
- Single-use



Qty: Various/Bx
Brand: Ultracell®

Description	SKU #
Classic Nasal Pack, 4.5 cm with String, 15 mm Thick, 20/Bx	10100-C
Classic Nasal Pack, 8 cm with String, 20/Bx	10110-C
Classic Nasal Pack, 10 cm with String 15 mm Thick, 20/Bx	10115-C
Anatomical Nasal Pack, 8 cm with String, 15 mm Thick, 10/Bx	10120-C

Epistaxis Products

The Ultracell® brand of epistaxis sponge provides highly absorbent treatment of nosebleeds.

- Fluid absorption causes the packing to expand producing a tamponading effect
- Available in a variety of sizes: 3 cm anterior to 10 cm posterior
- Ideal for use in emergency rooms, clinics, and doctor's offices
- Latex-free
- Single-use



Qty: Various/Bx
Brand: Ultracell®

Description	SKU #
Epistaxis Pack, 10 cm with String, 10/Bx	10140-C
Epistaxis Slimpack, 8 cm with String, 20/Bx	10146-C
Epistaxis Pack, 5.5 cm with String, 20/Bx	10150-C
Anterior Epistaxis Pack, 3 cm with String, 20/Bx	10160-C

Non-Ophthalmic

Sinus Packs

A complete line of Ultracell® brand sinus packs specifically designed for use following sinus surgeries.

- Ideal packing following FESS (Functional Endoscopic Sinus Surgery)
- Used for post-operative packing or dressing of sinus cavity
- Available in a range of lengths and thicknesses
- Patented layered design allows for easy manipulation
- Latex-free
- Single-use



Qty: 20/Bx
Brand: Ultracell®


Description	SKU #
Sinus Ultra Slimpack with String 2.5 cm x 0.6 cm x 2 cm	20202-C
Sinus Pack with String 3.5 cm x 1.2 cm x 1.2 cm	20210-C
Sinus Slimpack with String 3.5 cm x 0.9 cm x 1.2 cm	20211-C
Sinus Ultra Slimpack with String 3.5 cm x 0.6 cm x 1.2 cm	20212-C

Ear Wicks

The Ultracell® brand provides both ear wicks and ear packs. Much smaller, compressed form allows straight-forward insertion into swollen ear canal. Pediatric sizes available.

Ear Wicks

- Perfect for use with treatments for otitis externa
- Help maintain antibiotics in the ear canal
- Also available fenestrated
- Latex-free
- Single-use



Qty: Various/Bx
Brand: Ultracell®

Description	SKU #
Ear Wick 9 mm Diameter x 15 mm, 50/Bx	30300-C
Pediatric Ear Wick 7 mm Diameter x 12 mm, 50/Bx	30301-C
Fenestrated Ear Wick 9 mm Diameter x 15 mm, 50/Bx	30305-C
Fenestrated Pediatric Ear Wick 7 mm Diameter x 12 mm, 50/Bx	30306-C
Long Ear Wick 9 mm Diameter x 24 mm, 10/Bx	30310-C
Pediatric Long Ear Wick 7 mm Diameter x 20 mm, 10/Bx	30311-C
Pediatric Long Fenestrated Ear Wick 7 mm Diameter x 20 mm, 10/Bx	30316-C

FOR COMPLETE PRODUCT INSTRUCTIONS, PLEASE REFER TO THE BVI IFU OR PRODUCT INSERT

Questions: bvimedical.com/customer-support

BVI and all other trademarks (unless noted otherwise) are property of BVI. BVI ©2021



BVI
bvimedical.com

1594865-06-01